



Care For Your Prosthesis

- Regular cleaning once in 15 days.
- Wash the prosthesis between your hands with a non-abrasive, fragrance and moisturizer free soap and luke warm water.
- Polishing is recommended every 6 months.
- Prosthetic eye should be replaced approximately every 5 years.



SANKARA EYE HOSPITAL[®]
YOUR VISION OUR MISSION
BENGALURU

PROSTHETIC (ARTIFICIAL)
EYE CENTRE



Contact Us

+91 8123566262

080-69038900

ocularist.blr@sankaraeye.com



Second Floor , C Block
Sankara Eye Hospital , Kundalahalli
gate, Bangalore 560037

**Extensive experience
and expertise in
Artificial eye, serving as
a trusted and reliable
clinic for patients from
more than a decade.**



Our Services

Customised Prosthesis

Discover a world of confidence and comfort with our meticulously crafted customised ocular prosthesis, providing a seamless and natural look.

Ready Made Prostheses

Experience the convenience and immediate restoration of your appearance with readymade ocular prostheses.

Scleral Shells

Our clinic creates highly precise and customized scleral shells for a comfortable and realistic fit.

Hollow Prosthesis

Experience the lightweight comfort and lifelike aesthetics with our expertly designed hollow prostheses, tailored to your unique needs.

Modification/ Polishing of Prosthesis

Embrace perfection with our modification and polishing services, ensuring your prosthesis achieves the ideal fit and flawless finish.



Ocularistry Clinic

www.sankaraeye.com

Ocularistry @ Sankara is unique in our ability to provide both stock & customised prosthetic solutions to patients.

The centre is a training & research hub looking to enhance the care to those requiring these services across India.

Our patients have been from India, Srilanka, Bangladesh, Nepal, USA, Malaysia, Nigeria, Iran, Iraq, Yemen and many more.



Prosthetic eyes are a very common cosmetic treatment for those who have lost an eye/or both to create a balanced facial appearance.

**This helps to improve your physical appearance.
Lifts your self-esteem and confidence and visual aesthetic.**

How do we make it?

- Through a combination of advanced techniques, we create customised ocular prostheses that look like your natural eye.
- They are made using dimensionally stable, biocompatible material-PMMA.
- We take 2-3 days to prepare customised prosthesis.
- Three day procedure includes counseling, evaluation, measurement, hand painting the prosthesis, and fitting and polishing.
- The fitting of prosthesis is not at all painful.

Why our clinic is unique?

- Our team consists well experienced oculoplasty consultants and expert team of ocularists.
- Holistic approach to patient care.
- Exceptional patient satisfaction.
- Comprehensive support including counseling to address the psychological aspects of wearing a prosthetic eye.

